

## STANDARD DATA CONVERSION MAPPING

### System Source: Wave

**Extraction Process:** Local server-based data will be extracted by a Planet DDS Conversion Technician. Cloud data, images, and documents will be provided via SFTP by the legacy PMS IT Team.

Wave Conversions tend to be a different format for each conversion. We need enough time to verify the format and make changes before the Final. Two trials are required two months before the final to account for changes in the format.

Patient Demographics	Status
Title	✗
First Name	✓
Last Name	✓
Middle Initial	✗
Marital Status	✗
Nickname	✗
Birthdate	✓
Full Address	✓
Patient ID	✓
Social Security Number	✓
Assigned Location	✓
Patient Comments – Alerts	✓
Gender Information	✓
Patient Status	✓
Home Phone	✓
Cell Phone	✓
Work Phone	✓
Email	✓
Relationships – Responsible Parties, Siblings (Converted as Family Members), Professionals (Dr/Dentists), Other Referrals.	✓
Patient Recalls	✓

**KEY**

Data Field Converted: ✓

Data Field Not Converted, Not Available: ✗

Responsible Party	Status
First Name	✓
Last Name	✓
Full Address	✓
Social Security Number	✓
Gender Information	✓
Birthdate	✓
Email	✓
Home Phone	✓
Cell Phone	✓
Work Phone	✓
Marital Status	✗

Person Contact Info	Status
Phone Numbers	✓
Email Addresses	✓
Physical Addresses	✓

Practice/Appointments	Status
Schedule Views	✓
Location	✓
Appointments (Future)	✓
Appointments (Past)	✓

Ledger History *	Status
Transaction History	✓
Contracts	✓
Aging	✓
AR Reports Match?	✓

\* This is a balance-forward conversion. It retains the ledger history for viewing purposes only and cannot be queried for reporting. You must obtain all reports prior to conversion to Cloud 9 from your previous software.

**KEY**

Data Field Converted: ✓

Data Field Not Converted, Not Available: ✗

Treatment Card	Status
TC Columns – Uses Cloud 9 default.	✗
TC Entries	✓
TC Notes	✓
TC Plan – Treatment plan text is currently converted into treatment card plan notes (to the right of the Steps field in the Treatment Plan tab of the Treatment Card).	✓

Organization & Insurance	Status
Employers (for Responsible Parties)	✓
Employer Addresses	✗
Employer Contact Info (Phone)	✗
Insurance Companies	✓
Insurance Addresses	✓
Insurance Contact Info (Phone)	✓
Insurance Comments	✗
Insurance Claims (ADA Claims)	✗
Patient Insurance Policies	✓

Organization	Status
Employees (Orthodontists & Staff)	✓
Employee Addresses (Orthodontists & Staff)	✗
Employee Contact Info (Orthodontists & Staff)	✗
Employee To-Dos	✗

Questionnaires	Status
Questionnaires (Setup)	✗
Question Categories (Setup)	✗
Questions	✗
Answers	✗
Historical Patient Questionnaires	✗
Historical Patient Questionnaire Answers	✗

**KEY**

Data Field Converted: ✓

Data Field Not Converted, Not Available: ✗

Setup Types	Status
Appointment Types	✓
Patient Statuses	✓
Relationship Types	✓
Employee Types	✓
Contact Info Types	✓
Transaction Types	✓

Scanned Documents	Status
Scanned Documents	✓

**KEY**

Data Field Converted: ✓

Data Field Not Converted, Not Available: ✗

**NOTE:** Historically, we have converted these items from this software. However, every version is different and some offices place information in different fields. We cannot guarantee that every field will convert until we see your trial data.